



Recombinant Rat ADP-ribosylation factor-like protein 13A (Arl13a)

Product Code	CSB-EP718325RA-B
Abbreviation	Arl13a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AXT3
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MFRLLTACWS RQKATEGKQR NVTIIVIGLD NSGKSHLVAA FQRLIPSKMH SEMRPTQTTL LLDDYQVSVY DLTGDLKGRE KWPNYAEAH GLVFVVDSSD IARIQEVKII LTRLMFDKRV SGKPILILAN KQDKKDALLP CDIIEYLLLE RLVNETKSMC RVEPCSTIRN LPKRHQPIV IGLRWLLAAI GDRYDELCTR HQTPSMSNIS SSKNIRGCGE RCSSDSLSTR GVGVRNNHR QHFEKRQHLE HRQHVEQRHF EKRQHFESRQ HLIQRSLDAR PLKPILQKDG FRIRPKKNMS VTFALDVIME EGECSRKIGA PNISKPCNSQ CYDTKTPAPS ADANLFKARR PKKRIETWDT EEMFLEDPRG EAFRSCVTGH SSRSSRNTQS LYFTG
Source	E.coli
Target Names	Arl13a
Protein Names	Recommended name: ADP-ribosylation factor-like protein 13A
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.