



Recombinant Rat Something about silencing protein 10 (Utp3)

Product Code	CSB-MP718329RA
Abbreviation	Utp3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AXX4
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MVKRSRRRGA AQWATVRAKA GRTATDENED DLGSPSPSPGD SSYYQDQVDE FHEARSRAVL AKGWNEVESG EEDDDEEEEV LPLNIDDEND EDGESSEEEEE DGEDDDGGSS VQSEAEASVD PSLSWGQRKR LYYD TDY GSK SRGRQSQQEV EEEEEEEEE AQVIQRRLTQ ALQEDDFGVA WVEAFAPVP QVDEAETRVV KDLAKVSVKE KKLMLRKESP ELLELIEDLK VKLTEVKDEL EPLIQLVEKG VIPPGKGSQY LKTKYNLYLN YCANISFYLI LKARRVPAHG HPVIERLVTY RNLINKLSVV DQKLSSEIRH LLTAKEGAGK KDLNPKAKLT KTKPKSAKQT DVNADLTEEP EFDEAALEFY KEMEDRPELK RKKEENSAEE QALEEQNAKR AITYQIAKNR GLTPRRKKID RNPRVKHREK FRRAKIRRRG QVREVRREEQ RYSGELSGIR AGVKKSIKLLK
Source	Mammalian cell
Target Names	Utp3
Protein Names	Recommended name: Something about silencing protein 10 Alternative name(s): Charged amino acid-rich leucine zipper 1 Disrupter of silencing SAS10 UTP3 homolog
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.