



Recombinant Rat Coiled-coil domain-containing protein 91 (Ccdc91)

Product Code	CSB-YP718338RA
Abbreviation	Ccdc91
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AY97
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MDDDDFGGFE AAETFDGEQG GNQAVSPAVP WATFPAVSGV RLSPASPELI LDHDRSSPSS GHLRSDAVIS SPDDTRADSS LVNQTISKVQ LQQSAHATHLN IPLFPLGLTD ESNHGALALE DEPEGPGVHV SNSQLRQKIS SLETKLKASE EEKQRIKKDV ESLMEKHSV L EKDFLKEKEQ DAVSFQARYR ELQEKHKQEL EDMRKAGHEA LSIIVDEYKA LLQSSVKQQL DAIEKQYVSA IEKQAHRCCE LLHAQHQRLL EVLDTEKELL KEKIQEALTQ QSQEQKETLG KCLQEEMQKN KETLESVAVKL EKEAMKDVIT KAVEEERENL EKVHAEEREM WKTEHARDQE RVAEAIQAAV QEQQRMSQEA VKAAIAEEQR RSEKAMEEAV KRTRDELVEY VREQRRLDQV TRQRRSLSSLE LFLSCAQKQL SALIATEPVD IE
Source	Yeast
Target Names	Ccdc91
Protein Names	Recommended name: Coiled-coil domain-containing protein 91 Alternative name(s): GGA-binding partner
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.