



Recombinant *Xenopus laevis* FAS-associated factor 2-A (faf2-a)

Product Code	CSB-MP718352XBE
Abbreviation	faf2-a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AZH6
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MAALEERELS QEQTEKLLQF QDLTGIESMD QCRQTLQQHN WNIEAAVQDR LNEQEGVPRV FNNPPNRPLQ VNTADHRVYS YVVS RPQPRG LLGWGYLLIM LPFRITYYTL LDIFRFTLRF IRPDPRSRTV DPVGDVVSFI HLFEEKYGR I HPV FYQGTYS QALNDAKQEL RFLLVYLHGE DHQDSDDFCR NTLCTPEVTH FINSRMLFWA CSTNKPEGFR VSQALRENTY PFLGMIMLKD RRMTVVGRLE GLMQPQDLIN QLTFIIEANQ TYLVSERLER EERNETQVLR QQQDEAYLVS LRADQEKERK KKEKQEQKRR EEEEAQRKQM LEERKRNLE EEKERKSECL PAEPVPDHPD NVKIIFKMPN GTRVERRFLF TQSLSVIHDF LFSLKETPEK FQIVTSFPRR VLPCLPSEEI PVPPTLQEAG LSQSLLFVQ DLTDD
Source	Mammalian cell
Target Names	faf2-a
Protein Names	Recommended name: FAS-associated factor 2-A Alternative name(s): UBX domain-containing protein 8-A
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.