



Recombinant Debaryomyces hansenii Telomere replication protein EST3 (EST3)

Product Code	CSB-MP718421DIS
Abbreviation	EST3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BQG5
Product Type	Recombinant Protein
Immunogen Species	Debaryomyces hansenii (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (Torulaspora hansenii)
Purity	≥85% (SDS-PAGE)
Sequence	MNVNIKYPLV LHDNWLLES I TGQINKNVKY TNSLIKKDFI SPSSISKTP I LRILRFLKTT KSSDLTAILS DSTHTVIAIF PFDPAINFE IRYKHRITYH TPNSLILIKQ ANLRFVNNVE LTTEWGITIE DEIDVAVLEV LDLEIFQRDQ IMLGVNIENN LQLIYYDKSY LNLGKSRQR DSLGEKEIKE VISQSYDDVV SI
Source	Mammalian cell
Target Names	EST3
Protein Names	Recommended name: Telomere replication protein EST3
Expression Region	1-202
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.