



Recombinant *Debaryomyces hansenii* Ribosome biogenesis protein NSA1 (NSA1)

Product Code	CSB-MP718438DIS
Abbreviation	NSA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BUJ2
Product Type	Recombinant Protein
Immunogen Species	<i>Debaryomyces hansenii</i> (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (<i>Torulasporea hansenii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MKILVTSDDT GAVKEVCSR GTDTSKQDAK QPKSIKNICS EPGASVRTRV LHMVNFKSTF LVASRLGGSV TIYDLNHEEY ELVHTYQLAI SNEDKPISLI ELEEFDVVVV AFESGKVFVI NLNNGKFDLD PIEIQLPGKK GISAFINNP EAGVFACGGK ENDARIIRLF EGEITKDVFE AENKQEFFQP DVIFTARNVK NDHLDLRPPI WISSILFFEE KPKDGYKFLT STRYGQVRIY DTTHGKRPIQ DYKVCEKPIV TLNFADSEEE VIVSDTHNLV AKYSLAQIDS KAFKTHSASA GEITKPVSKL LGKFSAGGNT GAIVGVNIFD DEIISTGGLD RYLRTYDISS REILAKVYLG VQISDVLMLD SEDEEEEEPE SEMGDKRKKK IEKRRRKALE AEENESDEEE LWNQLEENNK TKIEKKKRR
Source	Mammalian cell
Target Names	NSA1
Protein Names	Recommended name: Ribosome biogenesis protein NSA1
Expression Region	1-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.