



Recombinant Debaryomyces hansenii Structure-specific endonuclease subunit SLX1 (SLX1)

Product Code	CSB-YP718448DIS
Abbreviation	SLX1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BWW4
Product Type	Recombinant Protein
Immunogen Species	Debaryomyces hansenii (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (Torulaspora hansenii)
Purity	≥85% (SDS-PAGE)
Sequence	MDDNEEYDTK RKTRCGLHVH PDFYGVYLLR SVPKPKSFYI GSTPNPQRRRL RQHNGELKNG GAYRTKKSGF RPWEMICLVY NFPSKNVALQ FEHALQHPYQ TRHIKSELRI THKRNSGNTL HHKLGNI RLL LGSSFFSRMG LKVLFLDPEV HSAWCINKFG VNVTDNVQVN VTRFEDYFSR DNNDSSGFS LSRQKESERI YFESSKKILF FDNQPCFICN ETIDYQSESE VSFSSKLDVD AYLREGNMPL LAICYHENC R KLYHLSCLGL HFLEKGDELS NKTDSEEKLG DTVNYLTPLQ GKCLSCNNFI NWAKLSKMST KLREYFLKDL LNTGATSQFI ENDADD
Source	Yeast
Target Names	SLX1
Protein Names	Recommended name: Structure-specific endonuclease subunit SLX1 EC=3.6.1.-
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.