



Recombinant *Yarrowia lipolytica* 40S ribosomal protein S6 (RPS6)

Product Code	CSB-YP718465YAP
Abbreviation	RPS6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6C169
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	>85% (SDS-PAGE)
Sequence	MKLNIAYP AQ GTQKCIDIED EHRV RGFY EK RMGQEVEADF LGDEFKGYVL KITGGNDKQG FPMKQGV MHP TRV RLLLSKD SSCYRSRRAG ERKRKSVRGC IVSSDLSVLS VVIVKQGDAD IEGLTTETVP RRLGPKRANN IRKLFALSKE DDVRQYVIRR EVTTKSGKTY TKAPKIQRLI TPRRLAHKAQ LREKKVKAIQ ASKDAAA EYA QLLAKRVAEK KSEKAEKKR RASSLR TQSV QA
Source	Yeast
Target Names	RPS6
Protein Names	Recommended name: 40S ribosomal protein S6
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.