



Recombinant *Kluyveromyces lactis* Chromatin structure-remodeling complex subunit SFH1 (SFH1)

Product Code	CSB-BP718582KBK
Abbreviation	SFH1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CSR7
Product Type	Recombinant Protein
Immunogen Species	<i>Kluyveromyces lactis</i> (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (<i>Candida sphaerica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MSLVPQAYLS NFHNRVKNE D VPLFTAAQLS RNKRGKQVNY AEFDTDLLDE FIDKNEDEDDL EDDVDDSDGR RRGGDYDQV DGSDGGAAAA AATAAAAAGL GEDGGIGEGT PGSGDPGAVT GGTPAADTGP DGTNDGTAGT SNPSSELDRI KLNDLPDLES QDDSATPLAL LKYPRIRETF VQSRIAISYK DLLGDSIQDP QQVDIESPIM VPIRLNVEFS GHKLADFFMW NLNDHSMTP E QFATILCQDL DFPVLSNPNN SPYTQIISMI NEQLQEYETL ASLQVPDLHV IINLTANLDS KLYDDTFEWN LNDDSLCPEQ FAELVVQDLG LQREFVPAIA HSLHESLLKV KKWLEGNLN LAHVENKAAF GYTSGIRLDI DTLGASWFPK VEVLSQWEIE KREIEKERNM RRLKRESAKV DDGRSRRRGK RRMDDLETTL RI
Source	Baculovirus
Target Names	SFH1
Protein Names	Recommended name: Chromatin structure-remodeling complex subunit SFH1 Alternative name(s): RSC complex subunit SFH1 SNF5 homolog 1
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.