



# Recombinant *Kluyveromyces lactis* ATP-dependent rRNA helicase RRP3 (RRP3)

<b>Product Code</b>	CSB-MP718586KKB
<b>Abbreviation</b>	RRP3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6CT85
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Kluyveromyces lactis</i> (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) ( <i>Candida sphaerica</i> )
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAKIEKKQAV KNNQVLSLAE KIKRKALEKQ QQAHANEPSP SDEDSAQSNS KDSNSNEQPE ESEEIFESFT ELDLPELIE ACKNLNYPKP TPIQSKAIPP ALKGSDIIGL AQTGSGKTAA FAIPILNQLW HDQQPYACI LAPTRELAQQ IKETFDSLGS LMGVIRSTCIV GGMSMMDQAR DLMRKPHEI ATPGRLMDHL ENTKGFNLRK LKYLVMDEAD RLLDMEFGPV LDRILNIPT QGRTTYLFSA TMTSKIDKLQ RASLTNPVKC AVSNKYQTVL TLVQTLMVVP GGLKNTFLIY LLNEFIGKST IVFTRTKANA ERISNLCNLL EFSATALHGD LNQNQRTGAL DLFKAGKRSI LVATDVAARG LDIPSVDIVI NYDIPVDSKS YIHRVGRTAR AGRSGKSISL VSQYDLELIL RIEDVLGKKL PKENVKDAI LALRDSVDKA NGEVVMELNR RNKEKQARGK GRRGRMATRD NMDREER
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RRP3
<b>Protein Names</b>	Recommended name: ATP-dependent rRNA helicase RRP3 EC= 3.6.4.-
<b>Expression Region</b>	1-487
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.