



Recombinant *Erwinia carotovora* subsp. atroseptica Pantothenate synthetase (panC)

Product Code	CSB-MP718646EAAB
Abbreviation	panC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6D1X5
Product Type	Recombinant Protein
Immunogen Species	Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) (<i>Erwinia carotovora</i> subsp. atroseptica)
Purity	≥85% (SDS-PAGE)
Sequence	MLIIETPLLL RREVRWRQE GKRIALVPTM GNLHDGHMTL VDEARARADI VIASVFNPM QFERPDDLAR YPRTLQEDCE KLNRRGADLV FAPSPDVIYP NGLESQTFVD VPGLSSMLEG ASRPGHFRGV ATIVSKLFNL VQPDLCFGE KDYQQLALIR QLARDMGYDI DIVGVPIVRA KDGLALSSRN GYLSAEERQL APQLYQLMMA LSAQLDNGDR QIDILLEQVA EKLREAGFTP DELFICDADT LQPLSAASTR AVILMAAWLG KARLIDNHQV DLTV
Source	Mammalian cell
Target Names	panC
Protein Names	Recommended name: Pantothenate synthetase Short name= PS EC= 6.3.2.1 Alternative name(s): Pantoate--beta-alanine ligase Pantoate-activating enzyme
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.