



Recombinant *Erwinia carotovora* subsp. atroseptica Acylphosphatase (acyP)

Product Code	CSB-EP718677EAAB
Abbreviation	acyP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6D6C4
Product Type	Recombinant Protein
Immunogen Species	Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) (<i>Erwinia carotovora</i> subsp. atroseptica)
Purity	≥85% (SDS-PAGE)
Sequence	MDKVATAAYV YGVVQGVGFR YSTQRQAAQL ELSGYARNCD DGSVEVVASG EPHAVEMLIE WLKQGGPRSA RVDKVLIEPH GKTDYKGFTI RY
Source	E.coli
Target Names	acyP
Protein Names	Recommended name: Acylphosphatase EC= 3.6.1.7 Alternative name(s): Acylphosphate phosphohydrolase
Expression Region	1-92
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.