



Recombinant Danio rerio Lipase member H (liph)

Product Code	CSB-YP718713DIL
Abbreviation	liph
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q6DBU8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	QECEEMT DLNFKDSLALG TSLKVRLLLY TRADPSCGQL LSHQEPFSNS QFNVSSVTTT LIHGYRPTGS PPVWMKQFVE FLLNRRDMNV IVVDWNRGAT NMNYWQVVKN TRKVANNLTD LIQKMKDNGA NLSSIHMIGV SLGAHISGFT GANFNGEIGR ITALDPAGPE FNGRPPEDRL DPSDALFVEA LHTDMDALGY RNLLGHIDYY ANGGADQPGC PKTILSGSEY FKCDHQRSVF LYMSSVNGSC PIIAYPCESY TDFQDGT CMD CGKFKSAGCP IFGYDSVRWR DTLVQLEQTR TYFQTNKASP FCKVGYKVDI VSWNQKTHWG YLTIKLSNGT EETQVELNHK SLKFERFQET SVLAQFERDI QPVKKITLKF CPRKGLRPRK KLRLHIRT PLQNHRLRPLC RYDLLLEESK DVTFKPIPCE DSNF
Source	Yeast
Target Names	liph
Protein Names	Recommended name: Lipase member H EC= 3.1.1.-
Expression Region	24-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>