



Recombinant Rat TANK-binding kinase 1-binding protein 1 (Tbkbp1)

Product Code	CSB-BP718746RA
Abbreviation	Tbkbp1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DG50
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MESMFEDDIS ILTQEALGPS EVWLDGPGDP SLGGDMCSAS HFALITAYGD IKERLGGLE ENATLRRRLK VYEIKYPLIT DFGEEHGFPL YEIKDGSLE VEKVSLQQLR NQFQHELQKN KEQEEQLGEM IQAYEKLCVE KSDLETELGE MRALVETHLR QICGLEKQLQ QQQGLRDAAF SSLSPPAVPA TACPDLDLHY LALRGGPALG HAGWPGPTSV SVSELERRL EEALEAAQGE ARGAQLREEQ LQAECERLQG ELKQLQETRA QDLASNQSEC GMAWVKRVGD DQVNLALAYT ELTEELGRLR ELSSLQGRIL RTLLQEARN AGQRHSPLSQ RHSPAPACPS PSPPARPPPC APCQSPAAQR RSPVPPCPSP QRRRSPASPS CSPVPQRRS PVPPSCQSPS PRRRSPVPPS CPAPQPRPPP PPGERTLAER AYAKPPSHA KAGFQGRSY SELAEGAAYA AASPAWLQAE AATLPKPRAY GGELYGPGRP LSPRRAFEGI RLRFEKQPSE EEWAMPASP PSPEASTIRC ASFCAGFPIP ESPAATAYAH AEHAQSWPSI NLLMETVGSD IRSCPLCQLG FPVGYRDDAL IKHIDSHLEN SKI
Source	Baculovirus
Target Names	Tbkbp1
Protein Names	Recommended name: TANK-binding kinase 1-binding protein 1 Short name= TBK1-binding protein 1 Alternative name(s): Protein ProSAPiP2
Expression Region	1-613
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.