



Recombinant Danio rerio UPF0705 protein C11orf49 homolog (zgc:92873)

Product Code	CSB-EP718749DIL-B
Abbreviation	zgc:92873
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DGK9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MKMNTERFNL TVDEYLAETN VLFYLNDAVT QLLEHKEEYT QFGVVRYFAE YFTSVKNGNH VLFREFSYIK ATPHNRESFI HIFWKCFRQI GKNGDLLAMS EYSSLLQLLC PDFPAEMVQN TARIVLIDDV TDCLMSFSDF IYSFQIQFFY EEFVESISAI YQDLLSGKNP NTVIVPTSTS VEQLSSSAND KPQVQEGVDA TIFCECIEGL CERFKHKYPS TVAIKEILDN TQRVSYFGFL MALAKHEGIN QDIGALPNKP DLLIDPEMDQ ELERLIAQVA ISPTSNNNSS SSALGQKEMS KKASPRKSLH QRKRIEMESD GSTEETDSSE N
Source	E.coli
Target Names	zgc:92873
Protein Names	Recommended name: UPF0705 protein C11orf49 homolog
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.