



Recombinant Danio rerio NF-kappa-B inhibitor-interacting Ras-like protein 1 (nkiras1)

Product Code	CSB-EP718750DIL
Abbreviation	nkiras1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DGQ1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MGKGCKVVVC GMASVGKTAI LEQLLYGSHT VGAETSDTQE DIYVASVETD RGVREQLRLY DTRGLREGLD LPKHFFSVAD GFVLVYSVDC LESFKKVEVL KKEIDRSRDK KEVMVMVLGN KCELRERRQV DQDTAQQWAR GEKVKLWEVT VTDRSTLIEP FTSLTSRLTQ PQSKSAFPLP GRKSKGTPSN DI
Source	E.coli
Target Names	nkiras1
Protein Names	Recommended name: NF-kappa-B inhibitor-interacting Ras-like protein 1 Alternative name(s): I-kappa-B-interacting Ras-like protein 1 Short name= Kappa B-Ras protein 1 Short name= KappaB-Ras1
Expression Region	1-192
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.