



Recombinant Rat Polycomb group RING finger protein 1 (Pcgf1)

Product Code	CSB-MP718777RA
Abbreviation	Pcgf1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DLV9
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MRLRNQLQSV YKMDPLRNEE EVRVKIKDLN EHIVCCLCAG YFVDATTITE CLHTFCKSCI VKYLQTSKYC PMCNIKIHET QPLLNLKLDL VMQDIVYKLV PGLQDSEEKR IRDFYQSRGL DRVSQPSGEE PALRGLGLPF TSHDHYRYD EQLSLCLERL SSGKDKNKNV LQNKYVRCSV RAEVRHLRRV LCHRLMLNPQ HVQLLDNEV LPDHMTMKQL WLSRWFGKPS PLLLQYSVKE KRR
Source	Mammalian cell
Target Names	Pcgf1
Protein Names	Recommended name: Polycomb group RING finger protein 1 Alternative name(s): Nervous system Polycomb-1 Short name= NSPc1
Expression Region	1-243
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.