



Recombinant Salmonella paratyphi A UDP-N-acetylmuramoylalanine--D-glutamate ligase (murD)

Product Code	CSB-EP718910SBM-B
Abbreviation	murD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PDC2
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MADYQDKNVV IIGLGLTGLS CVDFFLARGV TPRVMDTRVT PPGLDKLPQE VERHVGGGLND EWLLAADLIV ASPGIALAHP SLAAAASAGV EIVGDIELFC REAQAPIVAI TGSNGKSTVT TLVGEMAKAA GVNVGVGNI GLPALMLLDA DRELYVLELS SFQLETTSSL QAAAATVLNV TEDHMDRYPF GLQQYRAAKL RVYEKAKVCV VNADDALTMP VRGADERCVS FGVNMGDYHL NRQQGETWLR VKGEKVLNVK EMKLSGQHNY TNALAALALA DAVGLPRASS LKALTTFTGL AHRFQLALEH NGVRWINDSK ATNVGSTEAA LNGLHVDGTL HLLGGDGKS ADFSPLARYL TGDRIRLYCF GRDGAQLAAL RPEIAQQTET MEEAMRLLAP RVQPGDMVLL SPACASLDQF KNFEQRGDVF TRLAKELG
Source	E.coli
Target Names	murD
Protein Names	Recommended name: UDP-N-acetylmuramoylalanine--D-glutamate ligase EC=6.3.2.9 Alternative name(s): D-glutamic acid-adding enzyme UDP-N-acetylmuramoyl-L-alanyl-D-glutamate synthetase
Expression Region	1-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.