



Recombinant Salmonella paratyphi A Cyclic di-GMP regulator CdgR (cdgR)

Product Code	CSB-BP718939SBM
Abbreviation	cdgR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PH83
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MIASLDELYH SELFFLPVMD ENARLVGLEI IATFAAEDGA VRMPTELVP RLSVEEQYCL FVEKLALLET CQHFFIQHKL IAWLNLPPAI SDLLLLDSEL FSRAARFPFL ELAINENYPG LNQGKNNETL ANLAMHFPLM LANFGAGEAS TKAIFDGLFK RVMLDKNFIQ QRAEMISFEP FMHAIVAQIS SSCESLMIAG IDTEAMFARA APLGFSAFQG GLWPPVPVSQ LIKLVQR
Source	Baculovirus
Target Names	ydiV
Protein Names	Recommended name: Cyclic di-GMP regulator CdgR Short name= c-diGMP regulator
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.