



Recombinant Salmonella paratyphi A Histidine--tRNA ligase (hisS)

Product Code	CSB-MP719002SBM
Abbreviation	hisS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PNI1
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	AKNIQAIRG MNDYLPGETA IWQRIEGLTK NVLGSYGYSE IRLPIVEQTP LFKRAIGEVT DVVEKEMYTF EDRNGDSLTL RPEGTAGCVR AGIEHGLLYN QEQLRWYIGP MFRHERPQKG RYRQFHQLGA EVFGLQGPDI DAELIMLTAR WWRALGISEH VSLELNSIGS LEARANYRDA LVAFLEQHQE TLDEDCKRRM YTNPRLVLDS KNPDVQALLN DAPVLGDYLD DDSREHFTGL CKLLDAAGIA YTVNQRLVRG LDYYNRTVFE WVTNSLGSQG TVCAGGRYDG LVEQLGGRAT PAVGFAMGLE RLVLLVQAVN PEFIASPVVD IYLVAAGAQT QSAAMTLAER LRDEMPGVKL MTNHGGGNFK KQFARADKWG ARIALVLGES EVADGTVVVK DLRSGEQTAV AQDSVAAHLR TLLG
Source	Mammalian cell
Target Names	hisS
Protein Names	Recommended name: Histidine--tRNA ligase EC= 6.1.1.21 Alternative name(s): Histidyl-tRNA synthetase Short name= HisRS
Expression Region	2-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.