



Recombinant Rat Protein ZNF365 (Znf365)

Product Code	CSB-MP719019RA
Abbreviation	Znf365
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PQS2
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MQQTTFEESR YRWQDSLENV AVCLPFRCP R CGDHTRFRSL SSLRAHLEFS HSYEERTLLT KCSLLPSLKD TDLLRSSELP KQGKLLRGHA KVTQKPSYV NLYSISHGHS KDPKPFEMVA ERPVSYVQTY TAVDIRADSL DAPRSSGLP TQDTKAAFEA HVREKFNMV EAVDRRTIEKR IDKLTKELAQ KTAELLEVRA AFVQLTQKKQ EVQRRERALH KQVDVAVEMI AVLKQRLTES EEELLRKEEE VVTFNHFLEA AAEKEVQGKA RLQDFIENLL QRVELAEKQL EYYQSQQASG FSHDTSEHML TDISSSRKSR CLSRGHQHSV CNHPELKAHF HLGKRSYLKK AKDERAGMQP AKAIHEQAES PREFFRPAKK GEHLGLSRKG NFRPKMAKKK PTAIVNII
Source	Mammalian cell
Target Names	Znf365
Protein Names	Recommended name: Protein ZNF365
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.