



Recombinant Human DDB1- and CUL4-associated factor 10 (DCAF10)

Product Code	CSB-BP719034HU
Abbreviation	DCAF10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5QP82
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFPFGPHSPG GDGSAGAGAE EPTPHEGQAA ATGPPSPLHP GADATHPPPP ARSPRRPGAP SLSPAPRSGE LGLPGAPSS TASAPGEPSP PSPPCRRPGP DCRAKSRGRH GLGAGLGGPG ARLFGWLKER SLGRGLFVDP ARDNFRTMTS LYGSIHPADS VYLSTRTHGA VFNLEYS PDG SVLTVACEQT EVLLFDPISS KHIKTLSEAH EDCVNNIRFL DNRLFATCSD DTTIALWDLR KLNTKVCTLH GHTSWVKNIE YDTNTRLLVT SGFDGNVIIW DTNRYTEDGC PHKKFFHTRF LMRMRLTPDC SKMLISTSSG YLLILHDLDL TKSLEVGSYP ILRARRTTSS SDLTTSSSSS GPRVSGSPCH HSDSNSSEKH MSRASQREGV SPRNSLEVVT PEVLGESDHG NCITSLQLHP KGWATLLRCS SNSDDEECTC VYEFQEGAPV RPSRCSLR LTHYIEEANV GRGYIKELCF SPDGRMISSP HGYGIRLLGF DKQCSELVDC LPKEASPLRV IRSLYSHNDV VLTTKFSPTH CQIASGCLSG RVS LYQPKF
Source	Baculovirus
Target Names	DCAF10
Protein Names	Recommended name: DDB1- and CUL4-associated factor 10 Alternative name(s): WD repeat-containing protein 32
Expression Region	1-559
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.