



Recombinant Human Tetratricopeptide repeat protein 38 (TTC38)

Product Code	CSB-MP719094HU
Abbreviation	TTC38
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R3I4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AAASPLRDC QAWKDARLPL STTSNEACKL FDATLTQYVK WTNDKSLGGI EGCLSKLCAA DPTFVMGHAM ATGLVLIGTG SSVKLDKELD LAVKTMVEIS RTQPLTRREQ LHVSAVETFA NGNFPKACEL WEQILQDHPT DMLALKFSHD AYFYLGYYEQ MRDSVARIYP FWTPDIPLSS YVKGIIYSFGL METNFYDQAE KLAKEALSIN PTDAWSVHTV AHIHEMKAEI KDGLEFMQHS ETFWKDSDML ACHNYWHWAL YLIEKGEYEA ALTIYDTHIL PSLQANDAML DVVDSCSMY RLQMEGVSVG QRWQDVLPVA RKHSRDHILL FNDAHFLMAS LGAHDPQTTQ ELLTTLRDAS ESPGENCQHL LARDVGLPLC QALVEAEDGN PDRVLELLLP IRYRIVQLGG SNAQRDVFNQ LLIHAALNCT SSVHKNVARS LLMERDALKP NSPLTERLIR KAATVHLMQ
Source	Mammalian cell
Target Names	TTC38
Protein Names	Recommended name: Tetratricopeptide repeat protein 38 Short name= TPR repeat protein 38
Expression Region	2-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.