



Recombinant Pongo abelii Vacuolar protein sorting-associated protein VTA1 homolog (VTA1)

Product Code	CSB-EP719131PYX
Abbreviation	VTA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R5W5
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	AALAPLPPL PAQFKSIQHH LRTAQEHDKR DPVVAYYCRL YAMQTGMKID SKTPECRKFL SKLMDQLEAL KKQLGDNEAI TQEIVGCAHL ENYALKMFLY ADNEDRAGRF HKNMIKSFYT ASLLIDVITV FGELTDENVK HRKYARWKAT YIHNCLKNGE TPQAGPVGIE EDNDVEENED AGAASLPTQP TQPSSSSTYD PSNMPSGNYT GIQIPPGAHA PANTPAEVPH STGVASNTIQ PTPQTIPAID PALFNTISQG DIRLTPEDFA RAQKYCKYAG SALQYEDVST AVQNLQKALK LLTTGRE
Source	E.coli
Target Names	VTA1
Protein Names	Recommended name: Vacuolar protein sorting-associated protein VTA1 homolog
Expression Region	2-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.