



# Recombinant Pongo abelii Isocitrate dehydrogenase [NAD] subunit alpha, mitochondrial (IDH3A)

<b>Product Code</b>	CSB-MP719135PYX
<b>Abbreviation</b>	IDH3A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5R678
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	TGG VQTVTLIPGD GIGPEISAAV MKIFDAAKAP IQWEERNVTA IQGPGGKWMI PSEAKESMDK NKMGLKGPLK TPIAAGHPSM NLLLRKTFDL YANVRPCVSI EGYKTPYTDV NIVTIRENTE GEYSGIEHVI VDGVVQSIKL ITEGASKRIA EFAFEYARNN HRSNVTAVHK ANIMRMSDGL FLQKCREVAE NCKDIKFNEM YLDTVCLNMV QDPSQFDVLV MPNLYGDILS DLCAGLIGGL GVTPSGNIGA NGVAIFESVH GTAPDIAGKD MANPTALLLS AVMMLRHMGL FDHAARIEAA CFATIKDGKS LTKDLGGNAK CSDFTEEICR RVKDLD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	IDH3A
<b>Protein Names</b>	Recommended name: Isocitrate dehydrogenase [NAD] subunit alpha, mitochondrial EC= 1.1.1.41 Alternative name(s): Isocitric dehydrogenase subunit alpha NAD(+)-specific ICDH subunit alpha
<b>Expression Region</b>	28-366
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.