



# Recombinant Pongo abelii Multiple inositol polyphosphate phosphatase 1 (MINPP1)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP719174PYX                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Abbreviation</b>      | MINPP1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                 |
| <b>Uniprot No.</b>       | Q5R890                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen Species</b> | Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Sequence</b>          | SLLEPRDLVA SSSLSPYFGTK TRYEDVNPLL LSGPEAPWRD PELLEGSCPT<br>VQLVALIRHG TRYPTAKQIR KLRQLHGLLQ ARGSRDGGAG STGSRDLGAA<br>LADWPLWYAD WMDGQLVEKG RQDMRQLALR LASLFPALFS RENYGRRLRI<br>TSSKHRCMDS SAAFLQGLWQ HYHPGLPPPD VADMEFGPPT<br>VNDKLMRFFD HCEKFLTEVE KNATALYHVE AFKTGPEMQN ILKKVAATLQ<br>VPVNDLNADL LQVAFFTCSF DLAIKGVKSP WCDVFDIDDA KVLEYLNDLK<br>QYWKRGYGYT INSRSSCTLF QDIFRHLDKA VEQKQRSQPI SSPVILQFGH<br>AETLLPLLSL MGYFKDKEPL TAYNYKEQMH RKFRSGLIVP YASNLIFVLY<br>HCENAKTPKE QFRVQMLLNE KVLPLAYSQE TVSFYEDLRN HYKDILQSCQ<br>TSEECELARA NSTSDEL |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Target Names</b>      | MINPP1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Protein Names</b>     | Recommended name: Multiple inositol polyphosphate phosphatase 1 EC= 3.1.3.62 Alternative name(s): 2,3-bisphosphoglycerate 3-phosphatase Short name= 2,3-BPG phosphatase EC= 3.1.3.80 Inositol (1,3,4,5)-tetrakisphosphate 3-ph                                                                                                                                                                                                                                                                                                                     |
| <b>Expression Region</b> | 31-487                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                        |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.