



Recombinant Pongo abelii Protein TMED8 (TMED8)

Product Code	CSB-MP719201PYX
Abbreviation	TMED8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R9S0
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MSDLQAAEGP GSWSPARPG SAGGVGDCQG VEGSQAAASE NEDLENNKDT SLLASATDPE PCSPPHRPQM VSPGSKDATE DLRKATGALE AQALVKQDLL PADQAQVLNE MAKYQVPQRS GDIVMIQSEH TGAI DVLSAD LESADLLGDH RKVSPPLMAP PCIWTFKVK EFKSKLGKEK NSRLVVKRGE VVTIRVPTHP EGKRVCFEFA TDDYDIGFGV YFDWTPVTST DITVQVSDSC DDEDEEEEEEE EEIEEPVPAG DVERGSRSSL RGRYGEVMPV YRRDSDHRDVQ AGSHDYPGEG IYLLKFDNSY SLLRNKTLYF HIYYTS
Source	Mammalian cell
Target Names	TMED8
Protein Names	Recommended name: Protein TMED8
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.