



# Recombinant Pongo abelii Cytoplasmic dynein 1 light intermediate chain 2 (DYNC1LI2)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP719282PYX                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Abbreviation</b>      | DYNC1LI2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                        |
| <b>Uniprot No.</b>       | Q5RE09                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Immunogen Species</b> | Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Sequence</b>          | MAPVGVEKKL LLGPNGPAVA AAGDLTSEEE EGQSLWSSIL SEVSTRAGSK<br>LPSGKNILVF GEDGSGKTTL MTKLQGAEHG KKGRGLEYLY LSVHDEDRDD<br>HTRCNVWILD GDLYHKGLLK FAVSAESLPE TLVIFVADMS RPWTVMESLQ<br>KWASVLREHI DKMKIPPEKM RELERKFVKD FQDYMEPEEG CQGSPQRRGP<br>LTSGSDEENV ALPLGDNVLT HNLGIPVLVV CTKCDAVSVL EKEHDYRDEH<br>LDFIQSHLRR FCLQYGAALI YTSVKEEKNL DLLYKIYIVHK TYGFHFTTPA<br>LVVEKDAVFI PAGWDNEKKI AILHENFTTV KPEDAYEDFI VKPPVRKLVH<br>DKELAAEDEQ VFLMKQQSLL AKQPATPTRA SESPARGPSG SPRTQGRGGP<br>ASVPSSFPGT SVKKPDPNIK NNAASEGVLA SFFNSLLSKK TGSPGSPGAG<br>GVQSTAKKSG QKTVLSNVQE ELDRMTRKPD SMVTNSSTEN EA |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Target Names</b>      | DYNC1LI2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Protein Names</b>     | Recommended name: Cytoplasmic dynein 1 light intermediate chain 2<br>Alternative name(s): Dynein light intermediate chain 2, cytosolic                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Expression Region</b> | 1-492                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                               |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                                                                                                         |



of lyophilized form is 12 months at -20°C/-80°C.