



Recombinant Pongo abelii ATP synthase subunit b, mitochondrial (ATP5F1)

Product Code	CSB-BP719308PYX
Abbreviation	ATP5F1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RFH9
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	>85% (SDS-PAGE)
Sequence	PVPPLPEY GGKVRVYGLIP EEFFQFLYPK TGVTGPYVLG TGLILYALSK EIYVISAETF TALSLIGVMV YGIKKYGPVAV ADFADKLNEQ KLAQLEEAKQ TSIQQIQNAI DMEKSQQALV QKRHYLFDVQ RNNIAMALEV TYRERLYRVY KEVKNRLDYH IYVQNMMRQK EQEHMVNWVE KHVVQSISTQ QEKETIAKCI ADLKLLAKKA QAQPVM
Source	Baculovirus
Target Names	ATP5F1
Protein Names	Recommended name: ATP synthase subunit b, mitochondrial Short name= ATPase subunit b
Expression Region	43-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.