



Recombinant Pongo abelii Alpha-tocopherol transfer protein-like (TTPAL)

Product Code	CSB-EP719314PYX-B
Abbreviation	TTPAL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RFR0
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MSEESDSLRT SPSVASLSEN ELPPPPEPPG YVCSLTEDLV TKAREELQEK PEWRLRDVQA LRDMVRKEYP NLSTSLDDAF LLRFLRARKF DYDRALQLLV NYHSCRRSWP EVFNLNKPSA LKDVLASGFL TVLPHTDPRG CHVVCIRPDR WIPSNYPITE NIRAIYLTLE KLIQSEETQV NGIVILADYK GVSLSKASHF GPFIAKKVIG ILQDGFPIRI KAVHVVNEPR IFKGIFAIK PFLKEKIANR FFLHGSDLNS LHTNLPRSIL PKEYGGTAGE LDTATWNAVL LASEDDFVKE FCQPVPTCDS ILGQTLLEPEG LTSDAQCDLDS LRAVKSQLYS CY
Source	E.coli
Target Names	TTPAL
Protein Names	Recommended name: Alpha-tocopherol transfer protein-like
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.