



Recombinant Pig Carbonic anhydrase 3 (CA3)

Product Code	CSB-MP719334PI
Abbreviation	CA3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5S1S4
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	≥85% (SDS-PAGE)
Sequence	AKEWGYADH NGPDHWHELY PIAKGDNQSP IELHTKDIKH DPSLLPWTAS YDPGSAKTIL NNGKTCRVVF DDTYDRSMLR GGPLTAAYRL RQFHLHWGSS DDHGSEHTVD GVKYAAELHL VHWNSKYNSF ATALKHPDGV AVVGIFLKIG REKGEFQLVL DALDKIKTKG KEAPFTNFP SCLFPACRDY WTYHGSFTTP PCEECIVWLL LKEPITVSSD QMAKLRSLYS SAENEPPVPL VRNWRPPQPI KGRIVKASFK
Source	Mammalian cell
Target Names	CA3
Protein Names	Recommended name: Carbonic anhydrase 3 EC= 4.2.1.1 Alternative name(s): Carbonate dehydratase III Carbonic anhydrase III Short name= CA-III
Expression Region	2-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Carbonic anhydrase III (CAIII) is a member of a multigene family (at least six separate genes are known) that encodes carbonic anhydrase isozymes. These carbonic anhydrases are a class of metalloenzymes that catalyze the reversible hydration of carbon dioxide and are differentially expressed in a number of cell types. The expression of the CA3 gene is strictly tissue specific and present at high levels in skeletal muscle and much lower levels in cardiac and smooth muscle. A proportion of carriers of Duchenne muscle dystrophy have a higher CA3 level than normal. The gene spans 10.3 kb and contains seven exons and six introns.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.