



# Recombinant Human Tetratricopeptide repeat protein 39A (TTC39A)

<b>Product Code</b>	CSB-BP719385HU
<b>Abbreviation</b>	TTC39A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5SRH9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGQKGGHKDSL YPCGGTPESS LHEALDQCMT ALDLFLTQNF SEALSULKPR TKESMYHSLT YATILEMQAM MTFDPQDILL AGNMMKEAQM LCQRHRRKSS VTDSFSSLVN RPTLGQFTEE EIHAEVCYAE CLLQRAALTF LQGSSHGGAV RPRALHDP SH ACSCPPGPGR QHLFLLQDEN MVSFIKGGIK VRNSYQTYKE LDSLQSSQY CKGENHPHFE GGVKLGVGAF NLTL SMLPTR ILRLLEFVGF SGNKDYGLLQ LEEGASGHSF RSVLCVMLLL CYHTFLTFVL GTGNVNIEEA EKLLKPYLNR YPKGAIFLFF AGRIEVIKGN IDAAIRRFEE CCEAQQHWKQ FHHMCYWELM WCFTYKGQWK MSYFYADLLS KENCWSKATY IYMKAAYLSM FGKEDHKPFG DDEVELFRAV PGLK LKIAGK SLPTEKFAIR KSRRYFSSNP ISLPVPALEM MYIWNQYAVI GKQPKLTDGI LEITKAEM LEKGPENEYS VDDECLVKLL KGLCLKYLGR VQEAENFRS ISANEKKIKY DHYLIPNALL ELALLLMEQD RNEEAIKLE SAKQNYKNYS MESRTHFRIQ AATLQAKSSL ENSSRSMVSS VSL
<b>Source</b>	Baculovirus
<b>Target Names</b>	TTC39A
<b>Protein Names</b>	Recommended name: Tetratricopeptide repeat protein 39A Short name= TPR repeat protein 39A Alternative name(s): Differentially expressed in MCF7 with estradiol protein 6 Short name= DEME-6
<b>Expression Region</b>	1-613
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.