



Recombinant Human OTU domain-containing protein 3 (OTUD3)

Product Code	CSB-MP719399HU
Abbreviation	OTUD3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T2D3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSRKQAAKSR PGSGSRKAEA ERKRDERAAR RALAKERRNR PESGGGGGCE EEFVSFANQL QALGLKLREV PGDGNCLFRA LGDQLEGHSR NHLKHRQETV DYMIKQREDF EPFVEDDIPF EKHVASLAKP GTFAGNDAIV AFARNHQLNV VIHQLNAPLW QIRGTEKSSV RELHIAYRYG EHYDSVRRIN DNSEAPAHQ TDFQMLHQDE SNKREKIKTK GMDSEDDLDRD EVEDAVQKVC NATGCSDFNL IVQNLEAENY NIESAIIAVL RMNQGKRNN A EENLEPSGRV LKQCGPLWEE GSGGARIFGN QGLNEGRTEN NKAQASPSEE NKANKNQLAK VTNKQRREQQ WMEKKKRQEE RHRHKALESR GSHRDNNRSE AEANTQVTLV KTFAALNI
Source	Mammalian cell
Target Names	OTUD3
Protein Names	Recommended name: OTU domain-containing protein 3
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.