



# Recombinant Human MAGUK p55 subfamily member 7 (MPP7)

<b>Product Code</b>	CSB-MP719401HU
<b>Abbreviation</b>	MPP7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5T2T1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPALSTGSGS DTGLYELLAA LPAQLQPHVD SQEDLTFLWD MFGEKSLHSL VKIHEKLYHY EKQSPVPILH GAAALADDLA EELQNKPLNS EIRELLKLLS KPNVKALLSV HDTVAQKNYD PVLPPMPEDI DDEEDSVKII RLVKNREPLG ATIKKDEQTG AIIVARIMRG GAADRSGLIH VGDELREVNG IPVEDKRPEE IIQILAQSQG AITFKIIPGS KEETPSKEGK MFIKALFDYN PNEDKAIPCK EAGLSFKKGD ILQIMSQDDA TWWQAKHEAD ANPRAGLIPS KHQFERRLAL RRPEILVQPL KVSNRKSSGF RKSFRLSRKD KKTNKSMYEC KKSDQYDTAD VPTYEEVTPY RRQTNEKYRL VVLVGPVGVG LNELKRKLLI SDTQHYGVTV PHTTRARRSQ ESDGVEYIFI SKHLFETDVQ NKNFIEYGEY KNNYYGTSID SVRSVLAKNK VCLLDVQPHT VKHLRTLEFK PYVIFIKPPS IERLRETRKN AKIISSRDDQ GAAKPFTEED FQEMIKSAQI MESQYGHLFD KIIINDDLTV AFNELKTTTFD KLETETHWVP VSWLHS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	MPP7
<b>Protein Names</b>	Recommended name: MAGUK p55 subfamily member 7
<b>Expression Region</b>	1-576
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.