



Recombinant Human Glycosyltransferase 25 family member 3 (CERCAM)

Product Code	CSB-MP719403HU
Abbreviation	CERCAM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T4B2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AGVAESPL PAVVLAILAR NAEHSLPHYL GALERLDYPR ARMALWCATD HNVDNTTEML QEWLAAVGDD YAAVVWRPEG EPRFYPDEEG PKHWTKERHQ FLMEKQEAL TFARNWGADY ILFADTDNIL TNNQTLRLLM GQGLPVVAPM LDSQTYYSNF WCGITPQGYR RRTAEYFPTK NRQRRGCFRV PMVHSTFLAS LRAEGADQLA FYPPHPNYTW PFDDIIVFAY ACQAAGVSVH VCNEHRYGYM NVPVKSHQGL EDERVNFHIL ILEALVDGPR MQASAHVTRP SKRPSKIGFD EVFVLSLARR PDRRERMLAS LWEMEISGRV VDAVDGWMLN SSAIRNLGVD LLPGYQDPYS GRTLTKGEVG CFLSHYSIWE EVVARGLARV LVFEDDVRFE SNFRGRLERL MEDVEAEKLS WDLIYLGRKQ VNPEKETAVE GLPGLVVAGY SYWTLAYALR LAGARKLLAS QPLRRMLPVD EFLPIMFDQH PNEQYKAHFW PRDLVAFSAQ PLLAAPTHYA GDAEWLSDTE TSSPWDDDSG RLISWGSQSQK TLRSPRLDLT GSSGHSLQPQ PRDEL
Source	Mammalian cell
Target Names	CERCAM
Protein Names	Recommended name: Glycosyltransferase 25 family member 3 Alternative name(s): Cerebral endothelial cell adhesion molecule
Expression Region	23-595
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.