



# Recombinant Human Glycosyltransferase 25 family member 3 (CERCAM)

<b>Product Code</b>	CSB-YP719403HU
<b>Abbreviation</b>	CERCAM
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5T4B2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AGVAESPL PAVVLAILAR NAEHSLPHYL GALERLDYPR ARMALWCATD HNVDNTTEML QEWLAAVGDD YAAVVWRPEG EPRFYPDEEG PKHWTKERHQ FLMEKQEAL TFARNWGADY ILFADTDNIL TNNQTLRLLM GQGLPVVAPM LDSQTYYSNF WCGITPQGY Y RRTAEYFPTK NRQRRGCFRV PMVHSTFLAS LRAEGADQLA FYPPHPNYTW PFDDIIVFAY ACQAAGVSVH VCNEHRYGYM NVPVKSHQGL EDERVNFHIL ILEALVDGPR MQASAHVTRP SKRPSKIGFD EVFVISLARR PDRRERMLAS LWEMEISGRV VDAVDGWMLN SSAIRNLGVD LLPGYQDPYS GRTLTKGEVG CFLSHYSIWE EVVARGLARV LVFEDDVRFE SNFRGRLERL MEDVEAEKLS WDLIYLGRKQ VNPEKETAVE GLPGLVVAGY SYWTLAYALR LAGARKLLAS QPLRRMLPVD EFLPIMFDQH PNEQYKAHFW PRDLVAFSAQ PLLAAPTHYA GDAEWLSDTE TSSPWDDDSG RLISWSGSQK TLRSPRLDLT GSSGHSLQPQ PRDEL
<b>Source</b>	Yeast
<b>Target Names</b>	CERCAM
<b>Protein Names</b>	Recommended name: Glycosyltransferase 25 family member 3 Alternative name(s): Cerebral endothelial cell adhesion molecule
<b>Expression Region</b>	23-595
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.