



Recombinant Human Phosphatidylinositol 4-phosphate 5-kinase-like protein 1 (PIP5KL1)

| | |
|--------------------------|---|
| Product Code | CSB-MP719415HU |
| Abbreviation | PIP5KL1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q5T9C9 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MAAPSPGPRE VLAPSPEAGC RAVTSSRRGL LWRLRDKQSR LGLFEISPGH ELHGMTMMQ AGLWAATQVS MDHPPTGPPS RDDFSEVLTQ VHEGFELGTL AGPAFAWLRR SLGLAEEDYQ AALGPGGPYL QFLSTSKSKA SFFLSHDQRF FLKTQGRREV QALLAHLPRY VQHLQRHPHS LLARLLGVHS LRVDRGKKT Y FIVMQSVFYP AGRISERYDI KGCEVSRWVD PAPEGSPLVL VLKDLNFQGK TINLGPQRSW FLRQMELDTT FLRELVLDY SLLIAFQRLH EDERGPSSSL IFRTARVQVQ AQSPEESRAQ NRLLPDAPN ALHILDGPEQ RYFLGVVDLA TVYGLRKRLE HLWKTLRYPG RTFSTVSPAR YARRLCQWVE AHTE |
| Source | Mammalian cell |
| Target Names | PIP5KL1 |
| Protein Names | Recommended name: Phosphatidylinositol 4-phosphate 5-kinase-like protein 1 Short name= PI(4)P 5-kinase-like protein 1 Short name= PtdIns(4)P-5-kinase-like protein 1 EC= 2.7.1.68 |
| Expression Region | 1-394 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.