



Recombinant Human Phosphatidylinositol 4-phosphate 5-kinase-like protein 1 (PIP5KL1)

Product Code	CSB-YP719415HU
Abbreviation	PIP5KL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T9C9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAPSPGPRES VLAPSPEAGC RAVTSSRRGL LWRLRDKQSR LGLFEISPGH ELHGMTMMQ AGLWAATQVS MDHPPTGPPS RDDFSEVLTQ VHEGFELGTL AGPAFAWLR SLGLAEEDYQ AALGPGGPYL QFLSTSKSKA SFFLSHDQRF FLKTQGRREV QALLAHLPRY VQHLQRHPHS LLARLLGVHS LRVDRGKKT FIVMQSVFYP AGRISERYDI KGCEVSRWVD PAPEGSPLVL VLKDLNFQGK TINLGPQRSW FLRQMELDTT FLRELVLDY SLLIAFQRLH EDERGPSSS IFRTARVQV AQSPEESRAQ NRLLPDAPN ALHILDGPEQ RYFLGVVDLA TVYGLRKRLE HLWKTLRYPG RTFSTVSPAR YARRLCQWVE AHTE
Source	Yeast
Target Names	PIP5KL1
Protein Names	Recommended name: Phosphatidylinositol 4-phosphate 5-kinase-like protein 1 Short name= PI(4)P 5-kinase-like protein 1 Short name= PtdIns(4)P-5-kinase-like protein 1 EC= 2.7.1.68
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.