



Recombinant Human Glycolipid transfer protein domain-containing protein 1 (GLTPD1)

Product Code	CSB-MP719416HU
Abbreviation	GLTPD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TA50
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDDSETGFNL KVVLVSFKQC LDEKEEVLLD PYIASWKGLV RFLNSLGTIF SFISKDVVSK LRIMERLRGG PQSEHYRSLQ AMVAHELNSR LVDLERRSHH PESGCRTVLR LHRALHWLQL FLEGLRTSPE DARTSALCAD SYNASLAAYH PWVVRRAVTV AFCTLPTRV FLEAMNVGPP EQAVQMLGEA LPFIQRVYNV SQKLYAEHSL LDLP
Source	Mammalian cell
Target Names	CPTP
Protein Names	Recommended name: Glycolipid transfer protein domain-containing protein 1
Expression Region	1-214
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.