



# Recombinant Human U3 small nucleolar RNA-associated protein 14 homolog C (UTP14C)

<b>Product Code</b>	CSB-EP719418HU-B
<b>Abbreviation</b>	UTP14C
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5TAP6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNVNQVAENL ALSHQEELVD LPKNYPLSEN EDEGSDSDGER KHQKLEAII SLDGKNRRKL AERSEASLKV SEFSVSSEGS GEKGLGLADLL EPVKTSSSLA TVKKQLNRVK SKKVVELPLN KEKIEQIHRE VAFSKTSQVL SKWDPIILKN QQAEQLVFPL GKEQPAIPI EHALSGWKAR TPLEQEIFNL LHKNKQPVTD PLLTPMEKAS LQAMSLEEAK MHRAELQRAR ALQSYEAKA RKEKKIKSKK YHKVVKKGKA KKALKEFEQL QKVNPTVALE EMEKIENARM MERMSLKHQN SGKWAKSKAI MAKYDLEARQ AMQEQLAKNK ELTQKLQVAS EEEEEGGTE VEELLVPHVA NEVQMNVDGP NPWMFRSCTS DTKEAATQED PEQVPELAAH EVSASEAEER PVAEEEEILLR EFEERQSLRK RSELNQDAEP ASSQETKDSS SQEVLSELRA LSQKLKEKHQ SRKQKASSEG TVPQVQREEP APEEAEPDLL QRSERVQTL ELEEKGKEDC FQNKELPRPV LEGQQSERTP NNRPDAPKEK KEKEQLINLQ NFLTTQSPSV RSLAVPTIIE ELEDEEERDQ RQMIKEAFAG DDVIRDFLKE KREAVEASKP KDVDLTLPGW GEWGGVGLKP SAKKRRQFLI KAPEGPPRKD KNLPNVIISE KRNIHAAAHQ VQVLPYPFTH HRQFERTIQT PIGSTWNTQR AFQKLTPKV VTKPGHIIKP IKAEDVGYQS SSRSDLPVIQ RNPKRITTRH NKEEKL
<b>Source</b>	E.coli
<b>Target Names</b>	UTP14C
<b>Protein Names</b>	Recommended name: U3 small nucleolar RNA-associated protein 14 homolog C
<b>Expression Region</b>	1-766
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.