



# Recombinant Rat LIM domain-containing protein ajuba (Ajuba)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-BP719447RA   |
| <b>Abbreviation</b>      | Ajuba  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q5U2Z2   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Rattus norvegicus (Rat)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MERLGEKASR LLEKLRLSDS GSAKFGRRKG EASRSGSDGT PGAGKGRLSG<br>LGGPRKSGHR GANGGPGDEA LEPAREQGGL DAERNARGSF<br>EAQRFEGSFP GGPPPTRALP LPLSSPPDFR LETTAPALSP RSSFASSAS<br>DASKPSSPRG SLLLDGAGAS GAGGSRPCSN RTSGISMGYD QRHGSPLPAG<br>PCLFGLPLTT APSGYSSGGV PSAYPELHAA LDRLCAHRPV GFGCQESRHS<br>YPPALGSPGA LTGAVVGTAG PLERRGTQPG RHSVTGYGDC AAGARYQDEL<br>TALLRLTVAT GGREAGARGE PLGIEPSGLE ESPGSFVPEA SRSRIREPEA<br>REDYFGTCIK CNKGIYGQSN ACQALDSLYH TQCFVCCSCG RTLRCKAFYS<br>VNGSVYCEED YLFSGFQEAA EKCCVCGHLI LEKILQAMGK SYHPGCFRCI<br>VCNKCLDGVP FTVDFSNQVY CVTDYHKNYA PKCAACGQPI LPSEGCEDIV<br>RVISMDRDYH FECYHCEDCR MQLSDEEGCC CFPLDGHLLC<br>HGCHMQRLSA RQPPTNYI |
| <b>Source</b>            | Baculovirus  |
| <b>Target Names</b>      | Ajuba  |
| <b>Protein Names</b>     | Recommended name: LIM domain-containing protein ajuba  |
| <b>Expression Region</b> | 1-548  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.   |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.