



Recombinant *Xenopus tropicalis* Nucleoporin seh1 (seh1l)

Product Code	CSB-BP719461XBF
Abbreviation	seh1l
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U4Y8
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	MFVARSIAAD HKDLIHDVSF DFHGRRMATC SSDQSVKVWD KSENGDWHCT ASWKTHSGSV WRVTWAHPEF GQVLASCSFD RTAAVWEEIV GESNDKLRGQ SHWVKRTTLV DSRTSVTDVK FAPKHMGLML ATCSADGVVR IYEAPDVMNL SQWSLQHEIS CKLSCSCISW NPSSSRAHSP MIAVGSDDSS PNIMGKVQIY EYNENTRKYA KAETLMSVSD PVHDIAFAPN LGRSFHILAV ATKDVRIFTM KPLRKELSSS GGVTKFEIHT VAQFDNHNSQ VWRVSWNITG TVLASSGDDG TVRLWKANYM DNWKCIGVLK GDGNPVGNSC QGIFGSSVGS AIQSLQNSVN GTSSSSGRKHS
Source	Baculovirus
Target Names	seh1l
Protein Names	Recommended name: Nucleoporin seh1 Alternative name(s): Nup107-160 subcomplex subunit seh1
Expression Region	1-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.