



Recombinant Acanthamoeba polyphaga mimivirus Putative ankyrin repeat protein R760 (MIMI_R760), partial

Product Code	CSB-EP719509ADAZ-B
Abbreviation	MIMI_R760
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UPP7
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MDCFSQNEVF FMVTNTDEDI IIESSQKFFF TRIDNVDFDL KYGIYLREIT LPTNNTNFKL YTENGKLSAN IINFVVKHKL SDPETFYDLV EKGANIKPHL VNILFWASNH GFLNIIKYVH KYGIDLRTNN DYALRLACKG GHNNVVKFFI RNLDPKASN DFAICWAAEF GHYETVKILI KYNLDITTCS NYALRYACKN GYITLANYLI ENGADISAGD NYCIKIASEK GHYYMVNYLI SKGADV TAGN NYALNFAIKK DYIFIAKLLM KYGADIKEYL NDNLIDACRR NSIVTVRFLA ENGANINNSV VITTIVQHNY LDMLKCLVDT GLNLNSKKHD ILSIARDKGR TTIEYMMSQ NIDVLRWGIT NGYLDFVKMI PVIDNGFGIE LALILGHLDI AKYLISQNSN YQVDYDYVFR WSCRNGLTVI VKFLLSMKME PQSISYAVEW ATENGHLDIL KIILNKYMNV KIGKSIQTAR EKNYPEILDF LIGYLVKQEV D
Source	E.coli
Target Names	MIMI_R760
Protein Names	Recommended name: Putative ankyrin repeat protein R760
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.