



# Recombinant Acanthamoeba polyphaga mimivirus Putative ankyrin repeat protein R760 (MIMI\_R760), partial

<b>Product Code</b>	CSB-YP719509ADAZ
<b>Abbreviation</b>	MIMI_R760
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5UPP7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Acanthamoeba polyphaga mimivirus (APMV)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDCFSQNEVF FMVTNTDEDI IIESSQKFFF TRIDNVDFDL KYGIYLREIT LPTNNTNFKL YTENGKLSAN IINFVVKHKL SDPETFDYLV EKGANIKPHL VNILFWASNH GFLNIIKYVH KYGIDLRTNN DYALRLACKG GHNNVVKFFI RNLDPKASN DFAICWAAEF GHYETVKILI KYNLDITTCS NYALRYACKN GYITLANYLI ENGADISAGD NYCIKIASEK GHYYMVNYLI SKGADVTAGN NYALNFAIKK DYIFIAKLLM KYGADIKEYL NDNLIDACRR NSIVTVRFLA ENGANINNSV VITTIVQHNY LDMLKCLVDT GLNLNSKKHD ILSIARDKGR TTIEYMMSQ NIDVLRWGIT NGYLDFVKMI PVIDNGFGIE LALILGHLDI AKYLISQNSN YQVDYDYVFR WSCRNGLTVI VKFLLSMKME PQSISYAVEW ATENGHLDIL KIILNKYMNV KIGKSIQTAR EKNYPEILDF LIGYLVKQEV D
<b>Source</b>	Yeast
<b>Target Names</b>	MIMI_R760
<b>Protein Names</b>	Recommended name: Putative ankyrin repeat protein R760
<b>Expression Region</b>	1-501
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.