



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized protein L576 (MIMI_L576)

Product Code	CSB-EP719570ADAZ-B
Abbreviation	MIMI_L576
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UR96
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MDSSILSKLK YEPYSYSRCV KEYFITRDDN SDIKYFTSLI EIPTGNIGDN FLSFTDLNVF FYRLNSIDYY RSHKKFKICG AYITDLRSFK IKGVNDPTIN FNSFSFPQFT EFNPLIRSRN DRLFVVCVKFE CELSELYTLN AFMINFEGNY CEIPKDMYNC LNFSNQLEYF EGYIEMNYSE PIEYEFIVRA RSYLDIFCE FCLYDINFYM MNNDSEIVES LNDVAVLEME WENIDSRIDF KSFTEDNMLH LNTIEDVPMT NLPDVIPTDC FFSN
Source	E.coli
Target Names	MIMI_L576
Protein Names	Recommended name: Uncharacterized protein L576
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.