



Recombinant Human Forkhead box protein B2 (FOXB2)

Product Code	CSB-EP719630HU-B
Abbreviation	FOXB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VYV0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPRPGKSSYS DQKPPYSYIS LTAMAIQHSA EKMLPLSDIY KFIMERFPYY REHTQRWQNS LRHNLFSNDC FIKIPRRPDQ PGKGSFWALH PDCGDMFENG SFLRRRKRFK VLRADHTHLH AGSTKSAPGA GPGGHLHPHH HHHPHHHHHH HAAAHHHHHH HPPQPPPPP PPPPHMVHYF HQQPPTAPQP PPHLPSQPPQ QPPQQSQPQQ PSHPGKMQEA AAVAAAAAAA AAAAVGSVGR LSQFPPYGLG SAAAAAATAA ASTSGFKHPF AIENIIGRDY KGVQLAGGLP LASVMHHLGY PVPGQLGNVV SSVWPHVGVM DSVAAAAAAA AAAGVPVGPE YGAFGVPVKS LCHSASQSLP AMPVPIKPTP ALPPVSALQP GLTVPAASQQ PPAPSTVCSA AAASPVASLL EPTAP TSAES KGGSLHSVLV HS
Source	E.coli
Target Names	FOXB2
Protein Names	Recommended name: Forkhead box protein B2
Expression Region	1-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.