



Recombinant Human N-alpha-acetyltransferase 35, NatC auxiliary subunit (NAA35)

Product Code	CSB-BP719632HU
Abbreviation	NAA35
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VZE5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVMKASVDDD DSGWELSMPE KMEKSNTNWV DITQDFEEAC RELKLGELLH DKLFGLEFAM SAIEMMDPKM DAGMIGNQVN RKVLNFEQAI KDGTIKIKDL TLPELIGIMD TCFCLITWL EGHSLAQTVF TCLYIHNPDI IEDPAMKAFA LGILKICDIA REKVNKAAVF EEEDFQSMTY GFKMANSVTD LRVTGMLKDV EDDMQRRVKS TRSRQGEERD PEVELEHQQC LAVFSRVKFT RVLLTVLIAF TKKETSAAVE AQKLMVQAAD LLSAIHNSLH HGIQAQNDTT KGDHPIMMGF EPLVNQRLLP PTFPRYAKII KREEMVNYFA RLIDRIKTVK EVVNLNLHC ILDFFCEFSE QSPCVLSRSL LQTTFVLDNK KVFGLHLMQD MVKDALRSFV SPPVLSPKCY LYNNHQAKDC IDSFVTHCVR PFCSLIQIHG HNRARQRDKL GHILEEFATL QDEAEKVDAALHTMLLKQEP QRQHLACLGT WVLYHNLRLM IQYLLSGFEL ELYSMHEYYY IYWYLSEFLY AWLMSTLSRA DGSQMAEERI MEEQQKGRSS KKTKKKKKVR PLSREITMSQ AYQNMCGAMF KTMVAFDMDG KVRKPKFELD SEQVRYEHRF APFNSVMTTP PVHYLQFKEM SDLNKYSPPP QSPELYVAAS KHFFQAKMIL ENIPNPDHEV NRILKVAKPN FVVMKLLAGG HKKESKVPPE FDFSAHKYFP VVKLV
Source	Baculovirus
Target Names	NAA35
Protein Names	Recommended name: N-alpha-acetyltransferase 35, NatC auxiliary subunit Alternative name(s): Embryonic growth-associated protein homolog Protein MAK10 homolog
Expression Region	1-725
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.