



Recombinant Human Protein FAM22G (FAM22G)

Product Code	CSB-MP719633HU
Abbreviation	FAM22G
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VZR2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASNGAYPVL GPGVTVNPGT SLSVFTALPF ATPSPGPTHR PPLVTAVVPP AGPLVLSAFP STPLVAGQDG RGPSGAGASN VFVQMRTEVG PVKPPQAQTL ILTQAPLVWQ APGTLCGGVM CPPLLLLAAA PGVPVTSQAV VGGTQACEGG WSHGLPLPPP PPAAQVAPIV SPGNAGPWPQ GAHGEGLAP SQAKARPDDS CKPKSVYENF RLWQHYKPLA RRHLPQSPDT EALSCFLIPV LRSLARRKPT MTLEEGLWRA MREWQHTSNF DRMIFYEMAA KFLEFEAE MQIQKSQWMK GPQSLPPPAP PRLEPRGPPA PEVVKQPVYL PSKDGPKAPT ACLPPRPQR PAETKAHLPP PRPPRPAETK VPPEIPPEVV QEYVDIMEEL LGSHPGDTGE PEGQREKGV EQPQEEDGMT SDPGLLSYID KLCSQEDFVT KVEAVIHPRF LEELLSPDPQ MDLALSQEL EQEEGLTLAQ LVEKRLLSLK EKGCGRAAPR HGTARLDSSP SEFAAGQEA REVPDPQQRV SVETSPPQTA AQDPQGQGRV RTGMARSEDV AVLLGCQDSP RLKAVRPTSP PQDHRPTCPG LGTKDALGLP GESPVKESHG LAKGSSEETE LPGMVYVVG HHRLRPWRLS QSPVPSSGLL SPGGRGPQGA LQSPSAQKRG LSPSPSPASK SKKRPLFGSP SPAEKTPHPG PGLRVSGEQS LAWGLGGPSQ SQKRKGDPLA SRRKKKRHCS Q
Source	Mammalian cell
Target Names	NUTM2G
Protein Names	Recommended name: Protein FAM22G
Expression Region	1-741
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.