



Recombinant Legionella pneumophila Ribosomal RNA large subunit methyltransferase H (rlmH)

Product Code	CSB-YP719731LAAA
Abbreviation	rlmH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5WWX2
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila (strain Lens)
Purity	≥85% (SDS-PAGE)
Sequence	MLKITIITLG NKMPDWVNSG VNEYAKRFHD GIQIKLIEIP LLRRNKSSDL ARILEKESAL TKDALPANAR LIALDMLGKS FSSEELALKL TQLQQISSHL CFIIGGPEGL SNEILTLCD E RWSLSKLTLP HPLVRIILLE SLYRAWSIIN NHPYHK
Source	Yeast
Target Names	rlmH
Protein Names	Recommended name: Ribosomal RNA large subunit methyltransferase H EC= 2.1.1.177 Alternative name(s): 23S rRNA (pseudouridine1915-N3)-methyltransferase 23S rRNA m3Psi1915 methyltransferase rRNA (pseudouridine-N3-)-methyltransferase R
Expression Region	1-156
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.