



Recombinant Formimidoylglutamase (hutG)

Product Code	CSB-EP719804SMZ
Abbreviation	hutG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5X9K2
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MLEDYYPSTT SYYHGGIDDD LYTAKWGMVM TFLDLNDSSL TPFEGTHFAL IGFKSDKGVY INNGRVGAVE SPAAIRTQLA KFPWHLGNQV MVDVGNIDG PNRSLEQLQN SLSKAIKRMC DLNLKPIVLG GGHETAYGHY LGLRQSLSPS DDLAVINMDA HFDLRPYDQT GPNSGTGFRQ MFDDAVADKR LFKYFVLGIQ EHNNNLFLFD FVAKSKGIQF LTGQDIYQMG HQKICRAIDR FLEGQERVYL TIDMDCFSVG AAPGVSAIQS LGVDPNLAVL VLQHIAASGK LVGFDDVVEVS PPHDIDNHTA NLAATFIFYL VQIMAQHS
Source	E.coli
Target Names	hutG
Protein Names	Recommended name: Formimidoylglutamase EC= 3.5.3.8 Alternative name(s): Formiminoglutamase Formiminoglutamate hydrolase
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.